



CAMERA & CCTV ACCESSORY

<p>S Video Cable</p> 	<p>(#9901) S Video Cable</p> <p>Standard Length: 2m</p>
<p>BNC Video Cable</p> 	<p>(#9902) BNC Video Cable</p> <p>Standard Length: 2m</p>
<p>RCA Cable</p> 	<p>(#9903) RCA Cable</p> <p>Standard Length: 2m</p>
<p>S-Video to RCA adaptor</p> 	<p>(#9904) S Video to RCA adaptor</p>



BNC to RCA Adaptor



(#9905) BNC to RCA adaptor

USB Cable



(#9906) USB Interface Cable

Camera Power Adaptor



(#9907) Camera Power Adaptor 12V output

Camera Adaptor



(#9908) C/CS Camera Adaptor



VGA Cable



(#9909) VGA Cable

Video Splitter



(#9910) Video Splitter

Digital Camera Adaptor



(#9911) 10X Digital Camera Adapter Kit

- Suitable for 23.2mm and 30mm Eyepieces or C-Mount of microscope
- Lens with 28mm thread adapter for digital camera



(#9912) DSLR Digital Camera Adapter Kit



Digital Camera Adaptor



(#9913) 23.2/30mm Eyepiece Adaptor

(#9914) 23.2/30.5mm Eyepiece Adaptor

(#9915) DSLR Camera 2X Lens Adaptor

(#9916) Nikon T-Mount Adaptor

CCD Camera Adaptor



(#9917) 0.28x Video Coupler

(#9918) 0.45x Video Coupler

With standard C-Mount attached to CCD camera, the adaptor with 23.2mm tube, provides 30mm or 30.5mm adaptor to switch to your microscope

PRODUCT INFORMATION

For More Information

For more information on our products and services, kindly send an E-mail to sales@pinnacletech.biz

Tel: +(65) 6364 0626 Fax: +(65) 6894 9646