



# RING LIGHT ADAPTERS

## Product Information

## LED Ring Illuminator (ESD Version!)

TECHNICAL SPECIFICATION

	<p><b>Features.</b></p> <ul style="list-style-type: none"> <li>➤ ESD Version</li> <li>➤ Made of Black Anodized Aluminum</li> <li>➤ Designed &amp; Manufacture in Singapore</li> <li>➤ Different size available upon request</li> <li>➤ We provide customized Ring light Adapter upon your request</li> </ul>
--	--

Model Available	Type of Accessory Available*
PT-1501 Ring Adapter for Leica PT-1502 Ring Adapter for Olympus PT-1503 Ring Adapter for Nikon	<p>*Kindly Specify the type of model require when ordering</p>



## **PRODUCT INFORMATION**

Model	1501	1502	1503
Casing Material	Black Anodized, Al	Black Anodized, Al	Black Anodized, Al
Inner Diameter	Φ 66mm	Φ	Φ
Outer Diameter	Φ 59mm	Φ	Φ
Mounting Method	3 Thumb Screws	3 Thumb Screws	3 Thumb Screws

### **For More Information**

For more information on our products and services, kindly send an E-mail to [sales@pinnacletech.biz](mailto:sales@pinnacletech.biz)

Tel: +(65) 6364 0626 Fax: +(65) 6894 9646